

# Replacement Pull-String Parts

## Replacement Pull-String Parts

CS-01



✓ **Cord Stop -**

- Allows length of cord on Pull-String monitor to be adjusted
- For use with: TL-01, TL-02, TL-02E, TL-2000, TL-2200PS, TL-3000

GC-01



✓ **Garment Clip -**

- Allows cord to be clipped to the resident's clothing
- For use with: TL-01, TL-02, TL-02E, TL-02P, TL-2000, TL-2200PS, TL-3000

LC-02



✓ **Locking Garment Clip -**

- Allows cord to be clipped to the resident's clothing without slipping
- For use with: TL-01, TL-02, TL-02E, TL-02P, TL-2000, TL-2200PS, TL-3000

PS-01



✓ **Pull-String and Boot -**

- Replacement Pull-String and Boot
- For use with: TL-01

PS-02



✓ **Pull-String, Magnet & Clip -**

- Replacement Pull-String, Magnet & Clip
- For use with: TL-02, TL-02E, TL-2000, TL-2200PS, TL03000

PS-02P



✓ **Pull-String, Magnet & Clip -**

- Replacement Pull-String, Magnet & Clip
- For use with: TL-02P

PS-03



✓ **Snap Locking Garment Clip -**

- Allows cord to be clipped to the resident's clothing without slipping
- For use with: TL-01, TL-02, TL-02E, TL-2000, TL-2200PS, TL-3000

PS-04



✓ **Pull-String, Magnet & Clip -**

- Replacement Pull-String, Magnet & Clip
- For use with: TL-3000

RFC-01



✓ **Replacement Foam Cover for TL-2000 -**

- Protects casing from breaking and daily wear and tear
- For use with: TL-2000

RFC-02



✓ **Replacement Foam Cover for TL-3000 -**

- Protects casing from breaking and daily wear and tear
- For use with: TL-3000