

Low Cost Wireless Central Monitoring & Paging



433-MS
Motion
Sensor



Bed and
Chair Sensor
Pads



433-NCA
Wireless
Adapter for
Standard Monitors



Weight
Sensing
Floor Mats



433-EXT
Door/Window
Exit Sensor



433-NC
Nurse Call Button
(lanyard or wall mount)



433-CMU
Central Monitoring Unit



433-PGD
Caregiver Pager with
LCD Display



433-RB Caregiver Reset Button
(Allows caregivers to reset the
Central Monitoring Unit remotely)

Wireless
Signals

Wireless
Signal

Features:

- ✓ Monitors a large variety of wireless sensor devices at one location
- ✓ Numerical display and audible alert
- ✓ Suitable for residential care facilities and private homes

Also works with:



Caregiver Paging



Pager Case



AC Adapter for Wireless CMU



Pager Signal Bumper

Central Monitoring

Caregivers may easily monitor multiple residents at one location.

Caregiver Paging

Pagers allow caregivers to be aware of their resident's needs from anywhere within the facility.

Wireless Operation

No cords between devices reduces tripping and entanglement hazards, and allows facilities to reduce alarm noise in residents' rooms. No hard wiring required.

MFR#	SMART PRODUCT DESCRIPTION
433-CMU	Wireless Economy Central Monitoring Unit (1 year warranty)
433-PGD	Pager with LCD Display for use with Wireless CMU (6 month warranty)
433-NC	Nurse Call Button for use with Wireless CMU (90 day warranty)
433-RB	Wireless Reset Button for use with Wireless CMU (90 day warranty)
433-EXT	Door/Window Sensor for use with Wireless CMU (90 day warranty)
433-MS	Motion Sensor for use with Wireless CMU (90 day warranty)
433-NCA	Wireless Adapter for standard monitors (90 day warranty)
AC-04	AC Adapter for Wireless Central Monitor



SMART CAREGIVER CORPORATION

(800) 650-3637

www.smartcaregivercorp.com