



# Anti-Wandering Door System

## For a Facility of any Size!

The Smart DepartAlert™ Door Bar System functions with three easy-to-use components: a door strip monitor alarm, magnetic contact switch and resident wristband. When a resident wearing a wristband attempts to depart through the doorway, the door strip alarms audibly and visually.



### Install in 3 Easy Steps

- 1** Hang Monitor on Door
- 2** Plug in Monitor
- 3** Put Wristband on Resident



### Optional Central Monitoring



Mfr. #	Smart Product Description
TL-3005SYS	Single Door Bar Monitor System
TL-3004SYS	Double Door Bar Monitor System
TL-2012S	Resident Wristband Transmitter
TL-2014W	Wall Mount Reset Button
TL-2015REV-01	Central Monitor Unit
TL-2010WKP	Wireless Key Pad
DT-01	Door Tester & Wristband Control Unit
TL-2025	Magnetic Door Lock

