



Quiet, Cordless Fall Monitoring

Use Two Pads with One Monitor!



CordLess® Fall Monitor
PATENTED TECHNOLOGY

Wireless Signals

NEW!
Patented Technology

CordLess® Transmitter Inside Pad

CordLess® Transmitter Inside Pad

Features:

- ✓ Works with two sensor pads simultaneously
- ✓ Fully adjustable alert volume may be silent or any desired level with Chime Alert option
- ✓ Fully adjustable alert volume may be silent or any desired level with Chime Alert option
- ✓ **TamperProof™** technology (patented) with included Caregiver Key™
- ✓ Check Pad light alerts caregivers if a sensor pad's signal is lost
- ✓ Status Light flashes every 4 seconds when in use
- ✓ Nurse Call port, Low Battery light

No Cords Between Pads and Monitors

Reduce tripping hazards and eliminate broken or tangled pad cords.

Eliminate Bedside Alarms

Monitors may be placed outside resident's rooms and the alert may be adjusted to be quiet or any desired volume.

Works With Two Pads

Use one CordLess® monitor with two sensor pads for a safer and more complete resident monitoring environment.

MFR#	SMART PRODUCT DESCRIPTION
TL-2100GR1	CordLess® Fall Monitor (1 year warranty)
GBT-WI	CordLess® Bed Pad 20" x 30", (1 year pad life)
GBT-RI	CordLess® Bed Pad 10" x 30", (1 year pad life)
GCT-WI	CordLess® Chair Pad 10" x 15", (1 year pad life)
FMT-07C	CordLess® Floor Mat 24" x 48", (1 year warranty)



SMART CAREGIVER CORPORATION

(800) 650-3637

www.smartcaregivercorp.com