


Smart Unbreakable Magnet Pull-String Monitors with Unlimited Lifetime Warranty



- ✓ Reduce falls by alerting staff when a patient or resident attempt to get out of their chair.
- ✓ These Pull-String Alarm Monitors feature a cord that attaches via a garment clip to the residents clothing. When the resident attempts to get out of their chair, the pull-string disengages causing the fall alarm to sound, alerting the Caregiver of the resident's departure.
- ✓ Adjust the Cord Stop to desired Pull-String length for the resident's comfort and to prevent false alarms.



Unbreakable Magnet Pull-String Monitor Options



TL-2000

- ✓ **Unbreakable Magnet Pull-String Monitor**
 - Unlimited Lifetime Warranty
 - Uses a 9V Battery
 - Includes "unbreakable" cover for protection
 - Includes wheelchair mounting clip for easy mounting



TL-3000

- ✓ **Unbreakable Magnet Pull-String Monitor**
 - Unlimited Lifetime Warranty
 - Uses a 9V Battery
 - Includes "unbreakable" cover for protection
 - Adjustable volume and 3 tones
 - Alarm pause button silences an alarming monitor
 - Status light indicates when monitor is ready for use
 - Locking safety clip prevents residents from tampering
 - Safety release on cord separates when too much tension is on cord
 - Nurse call input
 - Includes wheelchair mounting clip for easy mounting

Mfr. #	Smart Product Description
TL-2000	"Unbreakable" Pull-String Monitor - Lifetime Warranty
TL-3000	"Unbreakable" Pull-String Monitor w/ Pause Button & Status Light - Lifetime Warranty

