


No more writing or checking dates on pads!



Monitor automatically indicates when to change pad



Change Pad Light

Change Pad Light

Indicates to caregiver when pressure pad needs replacement

Timer pauses when pressure is removed from the pad so only the days the pad is in use are counted

Tamper Proof

Switching to "Caregiver Key ON" mode on the side of the monitor means only the caregiver can silence the alarm

Reset/Status Button

Blinks when the system is actively working

SafeTCord™ Cord on Pads

Safety release cord prevents tripping hazards

Mfr. #	Product Description
TL-2100CP	Change Pad Indicator™ Monitor (patent pending)
TC-WI	Timed™ Chair Pad - 10" x 15" (1 year warranty) Coiled Cord
TC-90	Timed™ Chair Pad - 10" x 15" (90 day warranty) Coiled Cord
TC-45	Timed™ Chair Pad - 10" x 15" (45 day warranty) Coiled Cord
TB-WI	Timed™ Bed Pad - 20" x 30" (1 year warranty) Breakaway Cord
TB-RI	Timed™ Bed Pad 10" x 30" (1 year warranty) Breakaway or Straight Cord
TB-90W	Timed™ Bed Pad 20" x 30" (90 day warranty) Straight Cord
TB-90R	Timed™ Bed Pad 10" x 30" (90 day warranty) Straight Cord
TB-45W	Timed™ Bed Pad 20" x 30" (45 day warranty) Breakaway Cord

For use with
**Timed™ pressure
pads, floor mats,
and seat belts***

* The TL-2100CP monitor with Change Pad Indicator™ can be used with various Smart pressure pads, floor mats, and seat belts; however, the Change Pad Light will only function when the monitor is used with Smart Timed™ pressure pads.

Smart Caregiver Corporation • (800) 650-3637 • www.smartcaregivercorp.com

