



Economy Fall Monitoring



Chime Alert
Optional
Alert Tone

Features:

- ✓ Full range adjustable volume with optional **Chime Alert**
- ✓ **Pad Lost** technology informs caregivers if a pad becomes disconnected or stops working
- ✓ **On/Off** switch
- ✓ **Status Light** flashes every 4 seconds when in use
- ✓ **Nurse Call** port
- ✓ **Low Battery** light

Also works with:



Bed Sensor Pads



Floor Mats



Chair Belts

MFR#	SMART PRODUCT DESCRIPTION
TL-2100E	Economy Monitor (technology patent pending) (1 year warranty)
PPB-90W	Bed Sensor Pad 20"x 30", standard cord (90 day pad life)
PPC-45	Chair Sensor Pad 10"x 15", coiled cord (45 day pad life)
PPC-WI	Chair Sensor Pad 10"x 15", coiled cord (1 year pad life)

Low Cost

Economical and reliable mobility monitoring.

Universal Monitor

For bed pads, chair pads, floor mats and seat belts.

Status Light

When in use, the large status light blinks every 4 seconds.

Pad Lost Light

Technology informs caregivers if a pad becomes disconnected or stops working.

PATENTED TECHNOLOGY



SMART CAREGIVER CORPORATION

©2011 Smart Caregiver Corporation

(800) 650-3637

www.smartcaregivercorp.com