

**New!**

# No More Disturbing Alarm Noise!



**Musical Alert!**  
*Optional Alert Tone*



## Features:

- ✓ **Musical Alert** option with adjustable, full range volume control
- ✓ **Pad Lost** light alerts caregivers if a sensor pad is disconnected
- ✓ **Status Light** flashes every 4 seconds when in use
- ✓ **Nurse Call** port
- ✓ **Low Battery** light
- ✓ No "On/Off" switch

*Also works with:*



Chair Sensor Pads



Floor Mats



Chair Belts



Economy Central Monitor Unit

## Comes with Musical Alert

Residents and staff hear music rather than a startling alarm sound.

## Recordable Voice or Music

Easily customize the music or create a familiar voice recording.

## TamperProof™

In optional Caregiver Key Mode, the monitor can only be reset by the caregiver placing the Caregiver Key over the red dot – eliminating resident tampering.

**PATENTED TECHNOLOGY**



MFR#	SMART PRODUCT DESCRIPTION
TL-2100V	Melody Fall Monitor with Recordable Voice (1 year warranty)
PPB-RI, PPB-45	Bed Sensor Pad, 10"x 30", SafeTRelease™ Cord (1 year pad life)
PPB-90	Bed Sensor Pad, 5"x 30", Standard Cord (90 day pad life)
PPB-WI	Bed Sensor Pad, 20"x 30", SafeTRelease™ Cord (1 year pad life)
PPC-WI	Chair Sensor Pad, 10"x 15", Coiled Cord (1 year pad life)
FM-05, FM-07, FM-09	Floor Mat, 24"x 36", 24"x 48", 24"x 64", Gray (1 year warranty)
TL-2109	Easy Release Chair Belt (1 year warranty)
TI-2109V	Quick Release Hook & Loop Chair Belt (1 year warranty)



SMART CAREGIVER CORPORATION

(800) 650-3637

www.smartcaregivercorp.com