

# Pull-String Fall Monitors



## Features:

- ✓ DUAL Monitor works with pull-string or pressure pad
- ✓ Move easily from bed to wheelchair
- ✓ Adjustable alert volume
- ✓ Works with most existing nurse call light systems
- ✓ Non-Magnetic pull-string monitor available



*TL-3100 DUAL Works with Sensor Pads or as a Pull-String Monitor*

TL-3000  
Pull-String Monitor  
Pause Button,  
Auto-Reset



TL-02  
Pull-String  
Monitor



TL-02P  
Economy  
Pull-String  
Monitor



## Non-Magnetic Pull-String Monitor

Won't interfere with magnet-sensitive equipment



TL-01  
No Magnet  
Pull-String  
Monitor

MFR#	SMART PRODUCT DESCRIPTION
TL-3100	Dual Monitor: Pull String/Sensor Pad Monitor with nurse call port, status light, low battery, pad lost, adjustable tone & volume, (1 year warranty)
TL-3000	Pull-String Monitor w/Pause Button & Status Light, includes protective cover (10 year warranty)
TL-2200PS	Pull-String Monitor w/Pause Button & Status Light, (1 year warranty)
TL-02	Pull-String Monitor w/Magnet, (1 year warranty)
TL-02P	Pull-String Monitor w/Magnet, (6 month warranty)
TL-2000	"Unbreakable" Pull-String w/Magnet, (lifetime warranty)
TL-01	Non-Magnetic Pull-String Monitor, (1 year warranty)



SMART CAREGIVER CORPORATION

(800) 650-3637

www.smartcaregivercorp.com