



# Emergency Auto-Dialer



ED-01P

ED-01

## How It Works:



## Emergency Auto-Dialer Features:

- Opens Two-Way Communication Speakerphone
- No Monthly Fees
- Easy to program – plug into home telephone socket and use
- Supports up to 7 uniquely programmed help buttons
- AC adapter with battery back-up option
- Stores up to 4 different telephone numbers
- Can be used as a speaker phone
- Auto-redials until someone answers the phone call
- Phone line remains open after message is played allowing the caregiver to assess the situation
- Emergency Auto-Dialer has a One Year Warranty

Relatives can now have remote supervision of loved ones, without any monthly charges. This auto dialer opens up a two-way voice communication channel between the person at risk and the caregiver who receives the call. Dialer automatically dials up to 4 numbers of relatives, friends, neighbors or professional caregivers when the convenient wireless help button, worn as a pendant by the person being monitored, triggers the system. Once connected, a personal pre-recorded 20 second message giving caller information such as name, address and contact telephone number, is played. After, the channel remains open with a speakerphone allowing the caregiver to talk to the user and ascertain the nature of the caller's condition. The Emergency Auto Dialer may also be triggered from the main unit by pressing any one of the 3 orange "Emergency Call" buttons below the dial pad. The dialer buttons are large and suitable for sight impaired, and disabled people with limited hand movement.

Mfr. #	Product Description
ED-01	HelpSignal Emergency Auto Dialer
ED-01B	HelpSignal Emergency Auto Dialer Pendant



Demand Smart Products

**Smart Caregiver Corporation**  
 1229 N. McDowell Blvd., Petaluma, CA 94954  
 Toll Free: (800) 650-3637 • Fax: (707) 781-7440  
[www.smartcaregivercorp.com](http://www.smartcaregivercorp.com)

