



Pad Cords That Don't Break



SafeTRelease Cord™

- The SafeTRelease connector separates if the pad cord is under excessive tension
- Helps to reduce tripping, entanglement and choking hazards
- Durable, high quality cord construction
- Easily replaceable

Coiled Cords on Chair Pads

- Keeps chairs free from the clutter of loose cords
- Reduces potential tripping hazards
- Prevents cord damage

MFR. #	SMART PRODUCT DESCRIPTION
PPB-WI, PPB-90W, PPB- 45W	Bed Sensor Pad 20"x 30", SafeTRelease Cord - Warranties: 1 Year, 90 Day, 45 Day
PPB-RI, PPB-90, PPB-45	Bed Sensor Pad 10"x 30", SafeTRelease Cord - Warranties: 1 Year, 90 Day, 45 Day
PPC-WI, PPC-180, PPC-90	Chair Sensor Pad "10x15" - Warranties: 1 Year, 180 Day, 90 Day
TB-WI, TB-90W	"Timed" Bed Sensor Pad 20"x 30" SafeTRelease cord - Warranties: 1 Year, 90 Day
TB-RI, TB-90R	"Timed" Bed Sensor Pad 10"x 30", SafeTRelease Cord - Warranties: 1 Year, 90 Day
TC-WI, TC-90, TC-45	"Timed" Chair Sensor Pad 10"x 15" Coiled Cord - Warranties: 1 Year, 90 Day, 45 Day

Sensor pads work with all Smart TL-2100 Series Fall Monitors. "Timed" Pads work with all Smart TL-2100 Series Fall Monitors, but have a built in time tracking function that only works with the TL-2100CP Monitor.



SMART CAREGIVER CORPORATION

Tel: (800) 650-3637 www.smartcaregivercorp.com